



## MS Raw Data Conversion

### 1.1 HT-MS data extraction in Waters system by MassLynx


MassLynx software was used to extract data from Waters DI-QTOF MS acquisitions, generating data list files in .xls, .xlsx, or .csv formats containing  $m/z$  values and their corresponding intensities. The specific process is as follows:




① Open the TIC spectrum of the raw MS file



② Right click to select the range and extract the primary mass spectrum



③ Export raw data



④ Paste into an EXCEL spreadsheet

**Figure S1-1** HS-MS data extraction process using Masslynx

### 1.2 HT-MS data extraction in Agilent system by MassHunter

MassHunter software was utilized to extract MS data from Agilent DI-QTOF acquisitions, generating data list files in .xls, .xlsx, or .csv formats containing  $m/z$  values and their corresponding intensities. These files were subsequently imported into RapidMass for data processing, and the specific process is as follows:

### ③ Export MGF file




Figure S1-2 HS-MS data extraction process using Masshunter

## 1.3 MS data extraction in Thermo system and other systems by MZmine

MZmine software was employed to extract MS data from Thermo DI-Orbitrap acquisitions and DART-MS data from literature, producing data list files in .xls, .xlsx, or .csv formats containing  $m/z$  values and their corresponding intensities. Additionally, a "Multi-System Data Processing Guide" has been added to the user manual, which provides instructions on converting raw MS data into a universal format and offers solutions and recommendations for potential issues encountered during data processing across different systems. The complete user manual is available for download at <https://github.com/Katherine00689>.




Figure S1-3 HS-MS data extraction process using MZmine